

**screensafer**

# STAY SAFE

PERSONAL PROTECTIVE EQUIPMENT



**WEIGHT 55 GRAMS**



**ADJUSTABLE**

ELASTICS ALLOWS THE PERFECT FIT FOR ALL PEOPLE



**COMFORT**

WITH TPU MATERIAL REINFORCEMENT FOR COMFORTABLE FITTING



**CUSTOMIZABLE**


FRONT ALLOWS CUSTOMIZATION



**SAFETY AND HYGIENE**

DEVELOPED TO CE STANDARDS COMPLETELY CLOSED TOP TO PREVENT PARTICLES FROM ENTERING

REUSABLE



**CLEAR VIEW**

TRANSLUCENT MATERIAL ALLOWS A NATURAL VISION WITHOUT DISTORTIONS



**FRONT PROTECTION 150°**




**REINFORCEMENT**

TO KEEP ALWAYS THE BEST SHAPE AND BEST PROTECTION



**CLEAR PET FOIL**

EASY AND QUICK REPLACEMENT



screensafer

Developed in partnership with a team of doctors and medical specialists from hospitals and universities, we have designed a protection screen solution which minimize the risks of contamination.



We didn't reinvent the wheel but we have used our know-how and facilities to help in this crisis, and to accelerate the technical and productivity process in order to reach as fast as possible the help for the maximum number of people.


PRODUCT DEVELOPED BY:



[www.wdretail.pt](http://www.wdretail.pt)



[www.bsrp.eu](http://www.bsrp.eu)



screensafer

# STAY SAFE

PERSONAL PROTECTIVE EQUIPMENT



**STANDARD COLORS**

OTHER COLORS ON REQUEST



**PACKAGING**

10 UNITS PER BOX



## USAGE

The replacement of the **screensafer** components and the hygiene and disinfection procedures is the responsibility of the user. Always follow infection control procedures in accordance with the cleaning, disinfection and sterilization routines existing within each institution/company .

The **screensafer** must, essentially in the hospital environment and in contact with patients, be used in conjunction with other protective devices such as masks, caps and goggles and always in accordance with specific rules and instructions regarding any special procedure of infection control.


PRODUCT DEVELOPED BY:



[www.wdretail.pt](http://www.wdretail.pt)



[www.bspp.eu](http://www.bspp.eu)



screensafer

# STAY SAFE

PERSONAL PROTECTIVE EQUIPMENT



**STANDARD COLORS**

OTHER COLORS ON REQUEST



**PACKAGING**

10 UNITS PER BOX



## USAGE

The replacement of the **screensafer** components and the hygiene and disinfection procedures is the responsibility of the user. Always follow infection control procedures in accordance with the cleaning, disinfection and sterilization routines existing within each institution/company .

The **screensafer** must, essentially in the hospital environment and in contact with patients, be used in conjunction with other protective devices such as masks, caps and goggles and always in accordance with specific rules and instructions regarding any special procedure of infection control.

PRODUCT DEVELOPED BY:



[www.wdretail.pt](http://www.wdretail.pt)



[www.bspp.eu](http://www.bspp.eu)